

Product Bulletin

For Distribution To:

- International Distributors
- International Management
- U.S. Field Sales & Clinical Specialists
- U.S. Technical & Field Service

Date: October 1, 2013

RE: **Clean out procedure for switching from one type of disinfectant to another in the CER OPTIMA Reprocessor.**

NOTE: Please read through instructions thoroughly prior to performing this procedure. Verify that all the necessary supplies are at hand at the start of the process.

Remove the Existing Disinfectant

1. Dump existing disinfectant from the CER OPTIMA reservoir tank.
 - Press the HLD TO DRAIN key.
 - Press the START key.
2. Visually ensure that the reservoir tank is empty once the cycle is complete.
3. Remove the reservoir tank and wipe clean with a lint-free cloth.
4. Refit the tank in the cart or cabinet.
5. Remove the old disinfectant filter.
6. Install new disinfectant filter.

Rinse the CER OPTIMA with Water

1. Fill the reservoir tank with fresh water. This can be done either by directly pouring water into the reservoir, utilizing the transfer pump option or pouring the water into the basin and manually brining the water down into the reservoir. Refer to the CER OPTIMA User manual for instructions.
2. Run two HLD Rinse cycles with fresh water. NOTE: Water must be at the minimum required temperature for the selected HLD, otherwise the cycle will be halted due to low water temperature in the HLD phase of the cycle.
 - Select HLD RINSE Automatic Cycle and press START

After two cycles, dump the water from the reservoir tank.

- Press the HLD TO DRAIN Key
- Press the START Key.
- Remove the old disinfectant filter.
- Install new disinfectant filter.



MEDIVATORS Inc. Headquarters
14605 28th Avenue North
Minneapolis, MN
55447-4822 USA
Tel: +1.800.444.4729
Fax: +763.553.3387

MEDIVATORS BV
Sourthweg 11
6422PC Heerlen
The Netherlands
Tel: +31.45.5.471.471
Fax: +31.45.5.429.695

MEDIVATORS Asia/Pacific Pte Ltd
1 A International Business Park
#05-01 Singapore 609933
Tel: +65.6227.9698
Fax: +65.6225.6848

MEDIVATORS Beijing Representative Office
Room 708, 7th Floor Kaiheng Center, Block B
No. 2 Chaoyangmennei Street
Dongcheng District, Beijing China 100010
Tel: +8610.6567.8446
Fax: +8610.6567.8445

www.medivators.com

Product Bulletin

Rinse again the CER OPTIMA with Fresh water

1. Fill the reservoir tank with fresh water. This can be done either by directly pouring water into the reservoir, utilizing the transfer pump option or pouring the water into the basin and manually brining the water down into the reservoir. Refer to the CER OPTIMA User manual for instructions.
2. Run one HLD Rinse cycle with fresh water. NOTE: Water must be at the minimum required temperature for the selected HLD, otherwise the cycle will be halted due to low water temperature in the HLD phase of the cycle.
 - Select HLD RINSE and press START

After one cycle, dump the water from the reservoir tank.

- Press the HLD TO DRAIN Key
- Press the START Key.
- Remove the old disinfectant filter.
- Install new disinfectant filter.

Fill the CER OPTIMA with New Disinfectant solution

1. Fill the reservoir with new disinfectant solution. Refer to the CER OPTIMA User manual for information on reloading the disinfectant.
2. Verify the proper disinfectant profile is programmed



WARNING: When changing disinfectants, the proper disinfectant profile must be selected based on the type of disinfectant selected. Prior to patient use, always verify the CER OPTIMA is programmed correctly for the disinfectant being used. Call Medivators Technical Support at 800-444-4729 for details before proceeding.